



Beef Environmental Efficiency Programme - Suckler (BEEP-S)

DIRECTIONS

- 1) All faecal samples collected must be fresh.
- 2) Cows should be placed in a clean pen.
- 3) Using gloves, collect faecal samples from 10 different fresh faecal deposits and place them in 10 separate containers (You must collect at least a heaped teaspoonful of faeces).
- 4) Place all filled containers of faeces in the plastic zip lock bag.
- 5) Filled out submission form and make sure to include tag numbers.
- 6) Post samples, on the day of sampling, or the following day at most.
- 7) If samples are collected at the weekend store them in cool place preferable in a fridge not used for domestic purposes (DO NOT FREEZE or PLACE IN DIRECT SUNLIGHT) and post on Monday.
- 8) Cows should be dosed based on the result of the FEC result.
- 9) Applicants are advised to discuss the results of their FEC with their veterinary advisor.
- 10) We will forward on your details to the department.

KEY DATES

Eligible Animals: Beef
breed animals born
between 1 July 2019
and 30 June 2020

**CALL US ON (049) 854 1160 To arrange your
sampling kit.**

Testing, Kit & Postage all included for €25!!!!

Oldcastle Laboratories Ltd, Cogan St, Oldcastle, Co. Meath :
A82 HW90 : (049)8541160 : info@oldcastlelabs.ie